

*the new*



EMBRACE RF



MINIMALLY INVASIVE FACIAL  
CONTOURING (MINOR SURGICAL)

---

FACETite | ACCUTite | MORPHEUS8



# OFFER A FACIAL CONTOURING SOLUTION WHERE NO OTHER CURRENTLY EXISTS

EmbraceRF utilises InMode's innovative minimally invasive Radio Frequency Assisted Lipolysis (RFAL) technology to provide superior targeted lipolysis, delivering face sculpting results previously achieved only through more extensive surgical procedures.

The EmbraceRF is the perfect platform for aesthetic surgeons to offer a comprehensive, highly effective solution for gap patients who are not yet ready for a surgical procedure but want optimal outcomes, with minimal discomfort and downtime.

- » Subdermal adipose tissue remodelling of the face and neck
- » Sublative and ablative rejuvenation
- » In-office procedure that can be done under local anaesthesia
- » Non-excisional
- » Radio-Frequency Assisted Lipolysis (RFAL) technology to provide significant contraction
- » Real time measurements of skin temperature
- » Minimal downtime
- » Immediate results with continued improvements

## ACCUTite

## FACEtite



### NON - ASPIRATING

Cannula Length:	6 cm	10 cm
Cannula Diameter:	0.9 mm	1.3 mm
Thermal Sensor:	Dual	Dual/Single





EMBRACERF IS THE COMBINATION OF POWERFUL ENERGIES TO EFFECTIVELY TREAT THE SUBDERMAL ADIPOSE TISSUE OF THE FACE. THIS UNION RESULTS IN AN ENHANCED TREATMENT OUTCOME THAT CAN PROVIDE PHYSICIANS THE HIGHEST QUALITY FACIAL CONTOURING AVAILABLE.



## MORPHEUS8


Type:	Epidermal Resurfacing Sub-dermal Adipose Remodeling Device (SARD)
Configuration:	24 Pins
Treatment Depth:	Up to 4000+ $\mu$   4+ mm
Heating Depth:	Up to 5 mm
Output Pulse Power:	Up to 65 W
Frequency:	1 MHz
Repetition Rate:	Up to 1.5 pps

## IF A PATIENT HAS FAT AS WELL AS LOOSE SKIN, THE EMBRACE PROCEDURE CAN ADDRESS BOTH CONCERNS

EmbraceRF can remove the fat and contract the skin utilising InMode's proprietary **Facetite/Accutite** handpieces on the face or neck. Facetite/Accutite can be done under local anaesthesia in the office, with only two days of downtime. These non-aspirating cannulas allow physicians to work safely in the sub-dermal plane for tissue contraction and contouring. Utilising Radio-Frequency Assisted Lipolysis (RFAL) technology, these have been clinically proven in peer reviewed studies to provide significant contraction, with or without accompanying fat removal.

If the quality of the skin needs to be addressed, then external radio-frequency heat can be used to contour the face using **Morpheus8**, the unique RF Microneedling device for skin tightening, facial contouring and skin rejuvenation.

# INNOVATION THAT DELIVERS BEAUTY



Optional technology upgrades available

# BEFORE AND AFTER



PHOTOS COURTESY: DR M. CRAIG



PHOTOS COURTESY: DR J. BLOOM



PHOTOS COURTESY: DR P. ROVATTI



PHOTOS COURTESY: DR P. ROVATTI



AS A LEADER IN  
AESTHETICS , YOU CAN  
TRUST INMODE .

DELIVERING EXCELLENCE IN  
CLINICAL , TECHNICAL  
& MARKETING  
SUPPORT.



I N M O D E



InMode Australia



InMode\_Australia



InModeAustralia

INFO@INMODE.COM.AU | INMODE.COM.AU